

UNPACK SCOOTER AND RETAIN BOX FOR RETURN

- Lift Scooter out of box, swing main steering mast up to full extent. Push mast down onto locating peg, lock with quick release silver lever. **(Photo A)**
- Assemble the front wheel into the fork. Slide axle bolt through axle hole and tighten tri-knob.
- Fit right hand handle bar with brake lever onto main steering and lock into place. (Using black tri-knob) **(Photo B & C)**
- To set up for left or right leg use. **(Photo D)**
Move back axle away from the good leg of the user as indicated by the labels.
If this is not set for the correct non weight bearing leg – slacken both quick release levers on back axle and slide axle to the limit of the slot. Lock up quick release levers. Brakes will self-adjust.
- Adjust the knee post height so user's hips are level when non weight bearing leg is on the pad. Use the quick release silver lever on the knee post stem to do this. **(Photo E)**
- Adjust the knee pad sideways towards the user, using the two tri-knobs in the slots below the pad. Tighten tri-knobs. **(Photo F)**
- Adjust handle height to suit the user using clamp behind the basket hooks if needed.
- Brakes can be locked on by pushing the lever down away from the handgrip.
- Attach basket onto the hooks.
- To attach the Liquid Holster®, locate the mounting on the handlebar; attach the holder to this by inserting the small blue knob into the mounting. Quarter turn to lock.
- Attach crutches to crutch holder brackets, opposite side to user. Tips facing the front of Scooter, cuffs facing the rear. If back axle has been shifted, undo tri-knob on crutch holder. Swap crutch holder to opposite side then refit tri-knob and tighten. **(Photo G & H)**
- When using the Knee Scooter please exercise **CAUTION** when turning sharply to prevent tipping or overbalancing.
- DO NOT** overload the basket with heavy objects as this may cause the Knee Scooter to tip when turning.
- RETURNS** – Repack Knee Scooter into the box provided by reversing the above setup procedure.
See packaging list inside box and check that **ALL COMPONENTS** are packed.
- Contact Tubular Equipment Ltd, **Freephone 0800 757 588** so we can arrange freight pickup.
Please advise day for pickup and address if different from delivery address.



Photo A. Steering folding lock



Photo B. Handle bar with brake lever



Photo C. Handbrake lever on Scooter



Photo D. Back axle adjustment



Photo E. Knee post height adjustment



Photo F. Knee pad side adjustment



Photo G. Front crutch holder




Photo H. Back crutch holder